



All Trinity Fluids, LLC brands are designed to promote positive results, if used properly. You can choose to not deviate from our instructions or you can experiment with it based on your level of experience and desired personal look and feel you desire.

The following mixing instructions will give the embalmer a successful place to start.

AF Series mixing instructions:

- 120 – 180 mls of fluid per 4 liters of water.
- For trauma condition and difficult cases use 180 - 240 mls per 4 liters of water.
- All of our fluids are 6.8 – 7.2 pH. Giving a more natural feel. Review the firmness rating for your specific requirement.
- Complete fixation will occur 10 – 24 hours after embalming depending on cause of death, size of decedent and distribution of fluid.

Quick Reference Guide

Per 10 liters of water	50 kg	70 kg	90 kg	120 Kg
AF28	420 mls	450 mls	600 mls	750 mls
AF30 &	360 mls	420 mls	480 mls	600 mls
AF30 Rock'em	360 mls	420 mls	480 mls	600 mls
AF34	360 mls	420 mls	480 mls	600 mls
AMS Defiant	360 mls	420 mls	480 mls	660 mls